



Ocul-AR quick installation guide

1. Using your phone, navigate to www.bioimaging.fi/ocul-ar
2. Download the download Ocul-AR link (Requires Android 7.0 or newer)



3. During installation Your phone might ask if you are sure about installing:
 - since the source of the app is unknown, select yes.
 - If your phone does not allow installing, you have to navigate to settings and find the equivalent setting to let you install unconfirmed applications.
 - Your phone may ask if you want the application scanned and information sent, say no.
4. Ocul-AR should be now installed and accessible like any other application you have on your device.